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VORTEX Tracheo

Model: VORTEX Tracheo (Type 051)

VORTEX Tracheo

VORTEX Tracheoadapter

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VORTEX tracheostomy adapter

Read the instructions for use

Read these instructions for use carefully before using the product. Follow all instructions and safety instructions. Keep the instructions for use in a safe place.

Validity of instructions for use

Contact

E-mail: info@pari.de

Tel.: +49 (0)8151-279 220 (international)

+49 (0)8151-279 279 (German)

Copyright

©2022 PARI GmbH Spezialisten für effektive Inhalation
Moosstrasse 3, 82319 Starnberg, Germany. Technical and design changes and printing errors reserved. The same applies to illustrations.

Disclaimer

These instructions for use describe the components of PARI products and optional accessories. For this reason, these instructions for use also describe and illustrate features not present in your PARI product because they are, for instance, country-specific and/or optional. When using the systems, products and functions, the applicable country-specific regulations must be observed.

Trade marks

Registered trade marks of PARI GmbH Spezialisten für effektive Inhalation in Germany and/or other countries:

PARI®, VORTEX®

TABLE OF CONTENTS

IMPORTANT INFORMATION	5
Intended purpose	5
Indication	6
Contraindications.....	6
Labelling	6
Safety and warning instructions	7
PRODUCT DESCRIPTION	10
Components	10
Overview and product designations	10
Description of function	11
Material information	11
Operating life	12
USE	13
Function check	13
Preparing for treatment	13
Performing treatment.....	16
Ending the treatment	18
REPROCESSING IN A HOME ENVIRONMENT	19
Reprocessing cycles	19
Processing limits	19
Preparation	19
Cleaning	20
Disinfecting	22
Inspecting	23
Drying	24
Storage.....	24

REPROCESSING IN PROFESSIONAL HEALTH INSTITUTIONS.....	25
Reprocessing cycles	25
Processing limits	25
VORTEX with VORTEX tracheostomy adapter.....	26
FURTHER INFORMATION	31

1 IMPORTANT INFORMATION

1.1 Intended purpose

The VORTEX Tracheo is a holding chamber. It is used in conjunction with medication sprays, the so called metered dose inhalers for treatment of the lower airways.

The VORTEX Tracheo is equipped with a VORTEX tracheostomy adapter, with which treatment can be provided for tracheotomised patients in all age groups. The VORTEX tracheostomy adapter is designed only for use with a tracheostomy tube with a 15-mm-connector.

The VORTEX Tracheo must only be used after receiving professional instruction. The instruction can be performed by healthcare professionals or by someone who was instructed by healthcare professionals.

Individuals in the following groups must be supervised by a person who is responsible for their safety:

- Babies, infants, and children
- Individuals with limited capabilities (e.g. physical, mental, sensory)

If the patient is not able to operate the VORTEX Tracheo safely himself, the therapy must be carried out by the responsible person.

Only metered dose inhalers that have been approved for inhalation treatment must be used.

The VORTEX Tracheo can be used in a home environment and in professional health institutions. When used in a home environment, the VORTEX tracheo is intended for single-patient use only (no patient change). In a professional environment, the device can be used with different patients as long as the corresponding hygiene reprocessing measures are observed.

The VORTEX Tracheo is suitable only for patients who are able to breathe by themselves and are conscious.

The frequency and duration of use is determined by professional medical staff according to the individual needs of the patient.

1.2 Indication

For treatment of diseases of the lower airways in combination with a metered dose inhaler for use by tracheotomised patients.

1.3 Contraindications

There are no contraindications known to PARI GmbH.

1.4 Labelling

The following symbols can be found on the product and/or the packaging:

	Medical device
	Legal manufacturer
	Date of manufacture
	This product conforms to the EU Medical Device Regulation 2017/745.
	Consult instructions for use
	Item no.
	Production batch number, lot number
	Unique Device Identifier (UDI)
	VORTEX
	VORTEX tracheostomy adapter

1.5 Safety and warning instructions

The present instructions for use contain important information, safety instructions, and precautionary measures. The user must follow these in order to guarantee safe operation of this PARI product.

This PARI product must be used only as described in these instructions for use.

The instructions for use of the metered dose inhaler used must also be followed.

Labelling and classification of warning instructions

In these instructions for use, safety-critical warnings are categorised according to the following hazard levels:

DANGER

DANGER indicates a hazardous situation which will lead to very severe injuries or death if it is not avoided.

WARNING

WARNING indicates a hazardous situation which can lead to very severe injuries or death if it is not avoided.

CAUTION

CAUTION indicates a hazardous situation which can lead to mild or moderate injuries if it is not avoided.

NOTE

NOTE indicates a hazardous situation which can lead to material damage if it is not avoided.

General

If your health condition does not improve or it even worsens as a result of the treatment, seek professional medical advice.¹

Treatment of babies, infants, and anyone who requires assistance

For individuals who are not able to perform the therapy session without assistance or cannot appreciate the hazards, the risk of injury is greater. Such individuals include for example babies, children and people with limited capabilities. For these individuals, a person responsible for safety must supervise or carry out the application.

Hazard due to small parts which can be swallowed

The product contains small parts. Small parts can block the airways and lead to a choking hazard. Keep all components out of the reach of babies and infants at all times.

1) Professional medical staff: Doctors, pharmacists, and physiotherapists.

Hygiene

Comply with the following hygiene instructions:

- Use only components that have been thoroughly cleaned and dried. Contamination and residual moisture encourage bacterial growth, which increases the risk of infection.
- Before every use and reprocessing cycle, wash your hands thoroughly.
- Make absolutely sure you also carry out reprocessing before using the device for the first time.
- Always use drinking water for reprocessing.
- Make sure all components are dried properly after each reprocessing step.
- Do not keep the components in a damp environment or together with damp objects.

Reporting serious incidents


Report serious adverse incidents to the manufacturer and to the competent authority.

2 PRODUCT DESCRIPTION

2.1 Components

Please refer to the package for information on the supplied components.

2.2 Overview and product designations

(1)	Inhalation chamber (antistatic)
(1a)	Connecting ring for metered dose inhaler (not removable)
(2)	Mouthpiece with inspiratory and expiratory valve
(2a)	Protective cap
(3)	VORTEX tracheostomy adapter

2.3 Description of function

The VORTEX Tracheo helps the patient to inhale the medication from a metered dose inhaler evenly and completely.

The medication is sprayed into the inhalation chamber. There it is dispersed uniformly and can be inhaled completely by the patient in a few breaths. Consequently, it is not so important to coordinate the timing of releasing and inhaling the medication. It also prevents medication from being deposited in the wind-pipe. Side effects associated with premature deposition are avoided.

With the flexible connecting ring, it can be used with most common standard metered dose inhalers.

With the VORTEX tracheostomy adapter, the mouthpiece of the VORTEX Tracheo can be connected to a tracheostomy tube.

2.4 Material information

The individual components are made from the following materials:

Inhalation chamber	Aluminium
Mouthpiece	Polymethylpentene, polypropylene, thermoplastic elastomer
Inspiratory and expiratory valve in the mouthpiece	Silicone
Connecting ring for metered dose inhaler	Polypropylene, thermoplastic elastomer
Protective cap	Polypropylene
VORTEX tracheostomy adapter	Polypropylene

2.5 Operating life

The individual product components have the following expected lifetimes:

Product component	Service life
Inhalation chamber, mouthpiece (complete) VORTEX tracheostomy adapter	In a home environment[see: Processing limits, page 19]
Inhalation chamber, mouthpiece (complete) VORTEX tracheostomy adapter	In a professional environment[see: Processing limits, page 25]


3 USE

All the steps described below must be carried out correctly.

3.1 Function check

After receiving the VORTEX Tracheo, check that the path through the valve in the mouthpiece is unobstructed:

- Let the drinking water run out into the mouthpiece from the top.
The water should escape through the side openings in the mouthpiece.
- Let the drinking water run out through the bottom of the mouthpiece.
The water should run out of the lip attachment again.

- Dry the VORTEX Tracheo completely before using it for the first time.


If the valve is obstructed, the VORTEX Tracheo must not be used.

3.2 Preparing for treatment


Checking the valve position

Before every use, check that the valve is in the correct position inside the VORTEX mouthpiece:

- Check that the valve in the mouthpiece is in the position shown.

- If necessary, use a stream of water to correct its position.

- Dry the VORTEX Tracheo completely before using it.

Assembling the parts

WARNING

Choking hazard from breathing in small parts

The VORTEX Tracheo inhalation chamber and the VORTEX tracheostomy adapter are not completely closed. As a result, there is a danger that small parts can get into it, which can be breathed in during inhalation.

- Therefore, ensure that there are no foreign bodies inside the VORTEX Tracheo or the VORTEX tracheostomy adapter before every use.


CAUTION

Impaired treatment due to damaged or incorrectly assembled product components


Damaged components and an incorrectly assembled holding chamber may impair the function of the holding chamber and thus the treatment as well. Please also follow the instructions for assembly given below.

- Check all product components and the accessories before each use. Replace any broken, deformed, or seriously discoloured parts.


- Insert the mouthpiece in the inhalation chamber.

- Remove the cap on the metered dose inhaler.
- Shake the metered dose inhaler vigorously before every use.

- Insert the metered dose inhaler in the connecting ring on the inhalation chamber.

- Take the protective cap off the mouthpiece of the VORTEX Tracheo.
- Place the VORTEX tracheostomy adapter on the mouthpiece.

3.3 Performing treatment

All the safety instructions and warnings in these instructions for use must have been read and understood before any treatment is carried out.

CAUTION

Risk of underdosing with medication due to improperly assembled product components


If insufficient medication reaches the treatment site, the success of the treatment may be adversely affected.

- Make sure before the treatment that all product components are firmly connected to each other so that no medication can escape.


Unless otherwise prescribed by your doctor, the following numbers of breaths are enough to inhale the medication from the VORTEX:

- Adults need one breath
- Children older than 18 months need two to four breaths
- Children younger than 18 months need four to six breaths


Tip for visual assessment of breathing function by people helping patients to perform the therapy: The lateral valves must open visibly when exhaling.


- Assume a posture in which the upper body is upright (sitting or lying with the upper body raised).
- Connect the VORTEX tracheostomy adapter and the tracheostomy tube.


- Hold the VORTEX Tracheo as shown in the illustration.

- Actuate the metered dose inhaler at the start of a deep, slow intake of breath.
- Hold your breath briefly.

- Then breathe out slowly into the VORTEX Tracheo.


3.4 Ending the treatment

- Disconnect the VORTEX tracheostomy adapter together with the VORTEX from the tracheostomy tube.

⚠ CAUTION! Danger of injuring the tracheostoma!

Hold the tracheostomy tube firmly when removing the VORTEX tracheostomy adapter.

4 REPROCESSING IN A HOME ENVIRONMENT

4.1 Reprocessing cycles

VORTEX Tracheo (incl. VORTEX tracheostomy adapter)	<ul style="list-style-type: none">- Clean once a week. Clean once a day in the case of acute infection of the airways. Otherwise, whenever there is visible soiling.- Disinfect once per week. Clean once a day in the case of acute infection of the airways.
--	---


4.2 Processing limits

VORTEX Tracheo (incl. VORTEX tracheostomy adapter), disinfection	60 processing cycles, max. 1 year
--	-----------------------------------

4.3 Preparation

- Take the metered dose inhaler out of the connecting ring on the inhalation chamber and close it with the protective cap provided for this purpose.
- Take the VORTEX tracheostomy adapter off the mouthpiece of the VORTEX Tracheo.
- Disconnect the mouthpiece from the inhalation chamber.

Info: The protective cap must be open.

4.4 Cleaning

Precleaning

All individual parts must be pre-cleaned.

EQUIPMENT:

- Drinking water temperature of about 15 °C

PROCEDURE:

- Rinse all parts used for 2 minutes in running drinking water.


Manual cleaning

EQUIPMENT:

- Drinking water temperature of at least 40 °C
- Standard washing-up liquid²
- Receptacle with at least 3 l capacity

PROCEDURE:

- Add about 1 teaspoon of washing-up liquid to 3 l warm drinking water.
- Place all the parts in the washing-up water.

Application time: 5 minutes

- Occasionally move the parts back and forth.
- If there is visible soiling, use a medium-soft brush (such as a toothbrush) which is used exclusively for this purpose.

RINSING:

- Rinse all individual parts thoroughly under running drinking water at approx. 15 °C for 3 minutes.

DRYING:

- Shake the water out of all of the parts.
- Place all parts on a dry, clean and absorbent surface, and allow them to dry completely.

2) Validated with Palmolive®.

In the dishwasher

The individual parts can be cleaned in a standard household dishwasher, provided it is connected to a mains water supply of drinking water quality.

To ensure safety when handling the cleaning agent used, follow the corresponding instructions for use, particularly the accompanying safety instructions.

PROCEDURE:

i *Do not clean the individual components together with very dirty dishes.*

- Place all components in the crockery basket so that no water can collect in them.
- Select a program running at 50 °C minimum.

DRYING:

Ensure that there is no residual moisture remaining in the components. If necessary:

- Shake the water out of all of the parts.
- Place all parts on a dry, clean and absorbent surface, and allow them to dry completely.

4.5 Disinfecting

Disinfect all individual components after cleaning. Only components that have been cleaned can be disinfected effectively. The validated disinfection procedures are described below.

In boiling water

EQUIPMENT:

- Clean saucepan
- Drinking water

PROCEDURE:

CAUTION

Risk of infection due to moisture

Moisture encourages the growth of bacteria.

- Remove all parts from the saucepan as soon as the disinfection process is complete, and allow them to dry.

When disinfecting in boiling water, the mouthpiece must be detached from the inhalation chamber. Use enough water to ensure as little contact as possible between the mouthpiece and the inhalation chamber. The inhalation chamber becomes very hot and might deform the mouthpiece.

- Place all the individual parts in water at a rolling boil for at least 5 minutes. **NOTE! Risk of damage to plastic parts!** Plastic will melt if it comes into contact with the hot base of the pan. Make sure there is enough water in the pan to prevent the individual parts from touching the pan base.

DRYING:

- Shake the water out of all of the parts.
- Place all parts on a dry, clean and absorbent surface, and allow them to dry completely.

Using a standard thermal disinfecter for baby bottles (not a microwave oven)

EQUIPMENT:

- Thermal disinfecter with a runtime of at least 6 minutes

PROCEDURE:

CAUTION

Risk of infection due to inadequate disinfection

Inadequate disinfection encourages the growth of bacteria and thus increases the risk of infection.

- Make sure that the disinfecter is clean and operating properly before every disinfection process.
- Allow the disinfection to continue until the disinfecter switches off automatically or the minimum disinfection time stated in the instructions for use of the disinfecter has elapsed. Do not switch the device off prematurely.

CAUTION

Risk of infection due to moisture

Moisture encourages the growth of bacteria.

- Remove all parts from the disinfecter as soon as the disinfection process is finished, and allow them to dry.

Regarding the performance of the disinfection, the duration of the disinfection procedure and the quantity of water required, follow the instructions for use of the disinfecter you are using.

DRYING:

- After the disinfection process is complete, place all parts on a dry, clean and absorbent surface and allow them to dry completely. Or else leave all individual parts in the closed thermal disinfecter for max. 24 hours until the next use.

4.6 Inspecting

Inspect all product components after each cleaning and disinfection. Replace any broken, deformed or seriously discoloured parts.

4.7 Drying

After each cleaning and disinfection, place all components on a dry, clean and absorbent surface, and let them dry completely. The inhalation chamber can be dried with a clean, lint-free cloth.

4.8 Storage

Store the VORTEX Tracheo as described below:

- When it has dried completely, push the mouthpiece into the inhalation chamber.
- Place the protective cap over the mouthpiece.
- Store the product in a dry, dust-free place.

5 REPROCESSING IN PROFESSIONAL HEALTH INSTITUTIONS

5.1 Reprocessing cycles

Single patient use

VORTEX Tracheo (incl. VORTEX tracheostomy adapter)	<ul style="list-style-type: none">- Clean once a week. Clean once a day in the case of acute infection of the airways. Otherwise, whenever there is visible soiling.- Disinfect once per week. Clean once a day in the case of acute infection of the airways.
--	---

Before a change of patients

VORTEX Tracheo (incl. VORTEX tracheostomy adapter)	<ul style="list-style-type: none">- Cleaning- Disinfection- Sterilisation <p>Info: <i>Sterilisation is not necessary if a dedicated mouthpiece and VORTEX tracheostomy adapter are used exclusively by each patient.</i></p>
--	---

5.2 Processing limits

VORTEX Tracheo (incl. VORTEX tracheostomy adapter), sterilisation	30 processing cycles, max. 1 year
---	-----------------------------------

5.3 VORTEX with VORTEX tracheostomy adapter

Preparation

Dismantle the VORTEX Tracheo into its individual components.

Precleaning

All individual parts must be pre-cleaned.

EQUIPMENT:

- Drinking water temperature of about 15 °C

PROCEDURE:

- Rinse all parts used for 2 minutes in running drinking water.

Cleaning and disinfection

Generally, a mechanical process (instrument dishwasher) should be used for cleaning and disinfection.

NOTE

Damage to the product caused by cleaning agents and disinfectants that are incompatible with the material

Not all cleaning agents and disinfectants are suitable for use on the materials and product components. Unsuitable agents can damage the product and render it unusable.

- Before cleaning and disinfecting, check that the cleaning agent or disinfectant used is compatible with the materials of the product components that are to be processed.

Mechanical cleaning and disinfection

EQUIPMENT

- Cleaning and disinfection device in conformance with DIN EN ISO 15883
- Cleaning agent
- Neutralising agent if necessary

PROCEDURE

- Arrange the components so that they can be cleaned most effectively.
The inhalation chamber must be held upright.
- Select a suitable cleaning programme and agent.

VALIDATION

The method has been validated in Europe using:

- RDG G7836 and G7736 MCU produced by Miele (program Vario TD)
- Miele standard baskets (E142)
- Neodisher MediZym as a ph-neutral cleaning agent (0.5%, from Dr. Weigert)

Or:

- Neodisher Mediclean forte 0.5% (alkaline)
- For neutralising, use Neodisher Z (from Dr. Weigert)

DRYING

Even if the cleaner-disinfector is equipped with a drying function, make sure that no residual moisture is left in the components. If necessary, remove any residual moisture by shaking the parts and allow them to dry completely.

Chemical manual cleaning and disinfection

CAUTION

Risk of infection due to inadequate disinfection

Inadequate disinfection encourages the growth of bacteria and thus increases the risk of infection.

- Observe the specified mixing proportions and the specified treatment time.
- Ensure that all components are completely submerged in the solution for the whole of the treatment time. There must not be any air pockets or bubbles.

Alternatively, cleaning and disinfection can also be carried out manually. The validated agents are listed together with all relevant information in the following table:

	Korsolex® Endo Cleaner	Korsolex® Basic
Cleaning agent	×	–
Disinfectant	–	×
Contains aldehyde?	formaldehyde-free	×
Concentration	0.5%	4%
Duration	5 min.	30 min.
Temperature	50°C	max. room temperature

EQUIPMENT:

- Cleaning agent/Disinfectant
- Warm water
- Brush if applicable
- Basin

PROCEDURE:

If a disinfectant without a cleaning function is used, the following procedure must be performed twice: the first time with a cleaning agent, and then with the disinfectant.

To ensure safety when handling chemicals, follow the instructions for use of the agent used, particularly the accompanying safety instructions.

- Prepare the solution in the basin in accordance with the manufacturer's instructions. The quantity of solution required depends on the number of parts that are to be cleaned.
- Place all the individual parts in the prepared solution and leave them to soak.

Info: If the application period is exceeded significantly, the plastic parts may take on the smell of the disinfectant.

RINSING:

CAUTION

Risk of allergic reactions and irritation of the mucous membrane by disinfectants

Disinfectants can trigger allergic reactions or irritation of the mucous membrane on contact with the skin.

- Rinse the product thoroughly to ensure that no residues of the disinfectant remain on the PARI product.
- Rinse off all parts thoroughly in running drinking water.
Water temperature: approx. 15 °C, duration: 3 minutes

DRYING:

- Shake the water out of all of the parts.
- Place all parts on a dry, clean, and absorbent surface, and allow them to dry completely.

Sterilisation

CAUTION

Risk of infection by residual germs

If any germs are left on the parts that are to be sterilised, they can proliferate despite sterilisation. As a result, there is a danger of infection.

- Clean, disinfect, and dry all parts thoroughly before sterilising.
- Use only validated procedures for cleaning and disinfection.

EQUIPMENT

Steam steriliser – preferably with fractionated pre-vacuum – in accordance with DIN EN 285 or DIN EN 13060 (Type B)

PROCEDURE

- Pack all of the disassembled parts in a sterile barrier system in accordance with DIN EN 11607 (e.g., foil-paper packaging).

Maximum sterilisation temperature and holding time:

132 °C / 134 °C, at least 3 min.

VALIDATION

The method has been validated in Europe using: Steam sterilisers Euro-Selectomat and Varioklav 400E

Visual inspection and storage

Check all individual components. Replace any broken, deformed, or seriously discoloured parts.

- When it has dried completely, push the mouthpiece into the inhalation chamber.
- Place the protective cap over the mouthpiece.

Storage location:

- dry
- dust-free
- protected from sources of contamination

Optional: use sterile packaging

6 FURTHER INFORMATION

All product components may be disposed of with normal domestic waste. The country-specific disposal regulations must be observed.

All product components must be disposed of in accordance with the Ordinance on the European Waste Catalogue (EWC). Where applicable, country-specific disposal regulations must also be observed.

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©2022 PARI GmbH Spezialisten für effektive Inhalation, 051D4104 A it 2022-03-25