

Silicone Flat Drain

REF: 34.000.00.804 4mm for trocar CH 10
 34.000.00.807 7mm for trocar CH 15
 34.000.00.810 10mm for trocar CH 19

Unit: 5 pcs.

Packing Unit: 5 pcs.

Description

Silicone Flat Drain 4,7 and 10 mm with reservoir 100cm



Material Composition

Product: flat drain: silicone
 tube: silicone
 reservoir: silicone

Packing: individual packing: PE-Bag, paper, foil PA/PE
 case: cardboard
 outer case: corrugated cardboard

Sterilisation

EO-Gas