

Technical data

Gewicht Steuergerät	ca. 4,6 kg
Dimensions control unit	30 x 35 x 16 cm
Pressure power	30 - 240 mmHg
Time setting	5 - 99 min
Number of programs	8
Compatible with	Leg cuffs, Lymph-Flow Pants
Number of air chambers	12
Overlapping air chambers	✓
Deactivating air chambers	✓
Individually adjustable chamber pressure	✓
Quiet pump technology	65 db
HOLD-Timer	✓
INTERVAL-Timer	✓
Control	Steuergerät, Fernbedienung
SKU (Boots-Set)	VE-101
SKU (Lymph-Flow-Set)	VE-135
UVP (Boots-Set)	1.999 €
UVP (Lymph-Flow-Set)	2.299 €
Scope of delivery	Control unit, cuffs, hose set, operating instructions; power plug*.